

The VXR System's 2023 Redesign

It's here! The newly redesigned VXR (Veterinary X-Ray) System is available now, with customer feedback incorporated. The new system is retrofittable, compact, and highly practical for veterinary spaces.

Highlights

Compatible

We've removed the generator integration and added a "normal Bucky", so veterinarians can use the acquisition system of their choice.

Compact

The generator, available in 32kW with or without capacitor assist, has been relocated under the table.

Practical


A urine trap has been added, which is a common feature that techs appreciate.

Generator controls are conveniently located in front of the tube head and include an APR screen (pictured) for anatomical selection corresponding to generator techniques. These are pre-loaded with CR techniques, so *Fujifilm is providing a set of [DR APR settings](#) (available in the Box.com Service – VXR folder) that reduce the programmed dose by half.*

Techniques can also be customized and saved by the veterinary staff as needed.



Production units will have different accent color placement than pictured. No performance specifications are affected.



Configurations

The VXR System can be used with any detector system and is available with FDR SE Lite for customers on a small budget. Please contact your sales representative for pricing.

800051864	VXR System 32kW
800051865	VXR System 32kW Cap Asst
800052238	VXR 32kW System, FDR SE Lite C35
800052240	VXR 32kW System, FDR SE Lite C43
800052237	VXR 32kW Cap Asst, FDR SE Lite C35
800052239	VXR 32kW Cap Asst, FDR SE Lite C43

800052821	VXR Hand Switch
-----------	-----------------

Sales Resources

If needed, please request access via the [Fujifilm Distributor Resource Request Form](#).

The quote template, PowerPoint, and product photos are in the [Fujifilm Sales - VXR Veterinary X-Ray System](#) folder on Box.com.

The new system is installed and ready for customer visits in the Lexington, MA showroom.



Service Resources

If needed, please request access via the [Fujifilm Distributor Resource Request Form](#).

Service materials have been added to the [Fujifilm Service - VXR Veterinary X-Ray System](#) folder on Box.com, so anyone with existing access can also access the new materials.

The new system is straightforward to set up, and the manuals and SIBs are most likely sufficient for guiding personnel through the installation. As such, VXR system-specific service training has not been scheduled but can be if needed.

Discontinuations

This introduction officially marks the discontinuation of the original VXR and accessories.

800049887	VXR with FDR ES C35 - Integrated
800049888	VXR with FDR ES G35 - Integrated
800049104	VXR Veterinary X-Ray System - Integrated
800050006	VXR Free Cassette Wireless Kit
800049058	Compression Belt for VXR Table

Fujifilm appreciates your business and looks forward to providing imaging solutions to your organization. As always, please feel free to contact your channel sales support team with any questions or comments you may have.