



DIGITAL X-RAY

PODIATRY

SOLUTIONS


Revolutionary Imaging Systems




[2020imaging.net/podiatry](http://2020imaging.net/podiatry)

# Wireless and Tethered DR

 Cesium Scintillator

 16-bit Images



 Cassette Sized



 10x12"

## WIRELESS MOMENTUM DR PANEL

[2020imaging.net/momentum](http://2020imaging.net/momentum)





 Wireless Interface  
 Featherlight 3.75 lbs

 Lithium Ion Capacitor  
(7-yr warranty)  
 Rapid Charge

## TETHERED 1012V DR PANEL

[2020imaging.net/1012v](http://2020imaging.net/1012v)













 Tethered Interface  
 Lightweight 5.5 lbs

Detachable & Durable Cable

# New and Retrofit X-Ray

## NEW PXS710D X-RAY

-  High Frequency Output
-  Dose Efficient/Ext Tube Life
-  Superb Image Quality
-  Bi-directional Tube Movement
-  Integrated Cable Management
-  Enhanced Patient Support
-  UL Certified 450lbs weight limit
-  Qualifies for ADA Tax Credit
-  UL & CUL Listed
-  Health Canada Approved

New X-Cel X-Ray models also available

## RETROFIT ANY X-RAY

Wireless & Tethered DR Compatible  
All options qualify for ADA Tax Credit



\*example of X-cel MB-700



# Hospital-Grade Software

## ACQUISITION INTERFACE



- Under 6 second ⌚ Processing
- Background Submission
- Preset X-ray Exams
- Adult & Pediatric Settings
- Anatomical Program (APR)
- Fully Customizable
- Advanced Image Processing
- Image Enhancement Tools
- Touchscreen Enabled
- User Friendly Interface

## ADVANCED OPAL-RAD VIEWER



### Multi-fun-ton Studylist ⌵ Advanced Viewer ➡

- EHR Integration Ready\*\*
- Podiatric Toolset Available†
- Burn/Import CDs
- Fast, User Friendly Interface
- Custom Search/Query
- Image Fine Tuning
- Customized User Settings
- Exam Room Viewing



## MOBILE VIEWING WITH M-OPAL†

iPad®, iPhone®, Android™, Blackberry®, any mobile or smart device!



- Lightweight & Fast
- Patient List & Viewer
- Simple Annotations Enabled
- Tap, Drag, and Pinch
- Patient Filter & Search



## CLOUD-BACKUP INCLUDED\*\*

- Live Synchronization
- 448-bit Encryption (traffic)
- Enterprise Grade
- Fully Monitored
- Fully Customizable
- Recover Deleted Files



Advanced  
Superior  
Efficient  
Reliable  
Affordable

[2020imaging.net/podiatry](http://2020imaging.net/podiatry)

