

# ACCOAGENT PAUSCH


## VMS 100 Veterinary Elevating Table

Motorized elevating table top, with a travel range from approximately 8" to 48", allowing for easier movement of larger animal onto the table, and giving the operator easy access for positioning or any other point of care needs

Tabletop can move laterally  $\pm 5 \frac{1}{8}$ ", giving users a wider range of imaging area coverage for many different types of animals

The compact design is ideal for any facility size

The 3 -1 footswitch is tethered to the machine and can be operated from any place within range, giving users added convenience in their image positioning and care services



# ACCOAGENT PAUSCH

VMS 100



Veterinary Elevating Table