



site search

GO

HOME

ABOUT US

PARTS

PUBS

FINANCE

PRE-OWNED

MAIL-LIST

CONTACT US

SEARCH SITE

GO

NEW PRODUCTS

NEW! Zenith Embroidery ... Add your practice name, slogan, important wellness message to your upholstery in a beautiful and functional multi-colored embroidery.

All Zenith tables offer optional auto-cocking "AIR" drops.

SHOP BY CATEGORY

Chiropractic Tables

Automatic Cocking Drops

Elevation

Flexion Distraction

GONSTEAD

Hylo

Hylo w/ Elevation

Pediatric

Portable

Side Posture

Stationary

THOMPSON Tables

Upper Cervical

Decompression/Traction

Decompression

Intersegmental

Accessories

Education & Materials

Books, DVD, Tapes

Seminars

Electrotherapy & Ultrasound

Combination, Stim & Ultrasound

Laser & Light Therapy

Ultrasound Products

Iontophoresis Products

[Home](#) > [Chiropractic Tables](#) > [Hylo](#) > 220-Zenith Hylo Chiropractic Adjusting Table

220-Zenith Hylo Chiropractic Adjusting Table

[<< Previous in Hylo](#)[Next in Hylo >>](#)Price: **Call for Pricing****Order Today**Part Number: **220**[Add to a new shopping list](#)[Email this page to a friend](#)

This deluxe Zenith HYLO table incorporates both TERMINAL POINT cervical and pelvic drops to provide value with versatility for a wide range of techniques. You will feel the difference, and your patients will request the Zenith Hylo for their next treatment.

Features:

- 5 Year Limited Parts Warranty
- Tilt Headpiece
- Universal Headraiser raises or tilts
- Forward Motion Headpiece w/ lo-force drop mechanism
- Face Paper Toweling and Dispenser
- Soft Hinged Dorsal
- Swingaway Abdominal- Adjustable Tension
- Slide & Crank Pelvic
- Terminal Point Pelvic Drop w/ lo-force drop mechanism
- Adjustable Front Section
- Adjustable Ankle Rest Cushion
- Dual footswitches for easy access up, down, and stop
- Also front mounted up, down, an stop switches
- Heritage Upholstery in a variety of colors
- Base colors available in many coordinating colors
- Table Length 89"
- Table Height 22", 24", 26", or 28"
- Table Width 24"
- Cushion Width 18"
- Table Weight 390 lbs. depending on options

Additional Information[Click here ... UPHOLSTERY COLOR CHART](#)**Product Reviews**[Login to rate or review this product](#)**(0 Ratings, 0 Reviews)**

SPECIAL OFFERS

FEATURED PRODUCTS

Seminars-Dr. Rob Jackson




Posture Pump® Cervical Disc Hydrator® (Model 1000)



E44-Zenith Elevation-Series THOMPSON Table



Thompson Terminal Point Technique 8 DVD Set



Thumper Sport massager by THUMPER

Chronic Pain
Soft
Low-Volt Stim
T.E.N.S. Products

Laser & Light therapy
Massage Therapy
Equipment
Lotion

Physical Therapy Products
Braces & Supports
Therapeutic Exercise
Traction Therapy

Service Parts
Spine Tastic Support Belt
Supplies
Upholstery
Special Offers
Support Services

SHOPPING CART

Your cart is empty.



**BLOCKER SET - THOMPSON
TERMINAL POINT**

