

DJO GLOBAL

AUSTRALIA:

DJO Australia
PO Box 2057
Normanhurst
NSW 2076
AUSTRALIA
T: +1300 66 77 30
F: +1300 66 77 40
E: service@djortho.com.au

BENELUX:

DJO Benelux
Welvaartstraat 8
2200 Herentals
BELGIUM
T: Belgium 0800 18 246
T: Netherlands 0800 0229442
T: Luxemburg 8002 27 42
E: benelux.orders@DJOglobal.com

CANADA:

DJO Canada
6485 Kennedy Road
Mississauga
Ontario
L5T 2W4
CANADA
T: +1 866 866 5031
F: +1 866 866 5032
E: canada.orders@DJOglobal.com

CHINA:

DJO China
A312 SOHO ZhongShan Plaza
1055 W. ZhongShan Road
Shanghai, 200051
CHINA
T: (8621) 6031 9989
F: (8621) 6031 9709
E: information_china@DJOglobal.com

DENMARK, FINLAND, NORWAY & SWEDEN:

DJO Nordic AB
Murmansgatan 126
21225 Malmö
SWEDEN
T: Sweden 040 39 40 00
T: Norway 8006 1052
T: Finland 0800 114 582
T: Denmark +46 40 39 40 00
E: info.nordic@DJOglobal.com

FRANCE:

DJO France S.A.S.
Centre Européen de Fret
64990 Mouguerre
FRANCE
T: +33 (0)5 59 52 86 90
F: +33 (0)5 59 52 86 91
E: sce.cial@DJOglobal.com

GERMANY:

DJO Deutschland
ORMED GmbH
Merzhauser Str. 112
D-79100 Freiburg
GERMANY
T: +49 761 4566-01
F: +49 761 456655-01
E: infoservice@DJOglobal.com

ITALY:

DJO Italia Srl
Via Leonardo Da Vinci 97
Trezzaano Sul Naviglio
20090 Milano
ITALY
T: +39 02 484 63386
F: +39 02 484 09217
E: it.info@DJOglobal.com

SPAIN:

DJO Ibérica
Carretera de Cornellá 144, 1ª 4ª
Esplugues de Llobregat
08950 Barcelona
SPAIN
T: +34 943 638 167
F: +34 943 638 174
E: svc.cial@DJOglobal.com

SWITZERLAND:

Compex Médical SA
Ch. du Dévent
ZI Larges Pièces A
1024 Ecublens
SWITZERLAND
T: +41 (0) 21 695 2360
F: +41 (0) 21 695 2361
E: info@compex.ch

UK & IRELAND:

DJO UK Ltd
1a Guildford Business Park
Guildford
Surrey
GU2 8XG
UNITED KINGDOM
T: +44 (0)1483 459 659
F: +44 (0)1483 459 470
E: ukorders@DJOglobal.com

UNITED STATES:

DJO Global, Inc
1430 Decision Street
Vista
CA 92081-8553
U.S.A.
T: 1 800 336 6569
F: 1 800 936 6569

DJO GLOBAL, EXPORT CENTERS

ASIA-PACIFIC:

DJO Asia-Pacific Limited
Unit 1905, 19/F, Tower II
Grand Central Plaza
138 Shatin Rural Committee Road
Shatin
HONG KONG
T: +852 3105 2237
F: +852 3105 1444
E: info.asia@DJOglobal.com

EUROPE, MIDDLE EAST & AFRICA:

DJO Benelux
Welvaartstraat 8
2200 Herentals
BELGIUM
T: +32 (0) 14248350
F: +32 (0) 14248358
E: info.emea@DJOglobal.com

LATIN AMERICA:

DJO Global, Inc
1430 Decision Street
Vista
CA 92081-8553
U.S.A.
T: 1 800 336 6569
F: 1 800 936 6569
E: info.latam@DJOglobal.com



Together in Motion™

 **chattanooga™**

Intelect® Vet | 2013



INTRODUCTION



Chattanooga is the world's largest manufacturer of rehabilitation equipment for treating musculoskeletal, neurological and soft tissue disorders. For over six decades, Chattanooga has set a benchmark for leadership, reliability and excellence. Our products contribute to better treatment outcomes in hospitals, clinics and home settings worldwide. We lead the physiotherapy industry by example and, through continuous innovation, provide real world solutions for clinicians and their patients.

As the world's leading brand in rehabilitation, Chattanooga is sold in over 80 countries worldwide. Chattanooga is directly represented in our key markets but also has long lasting relationships with dedicated distributors in the USA, Canada, Latin America, EMEA and Asia Pacific. Our Chattanooga distributors and customers are supported for order placement, training, product delivery and after sales service by a dedicated team of highly skilled DJO employees based all over the world.

Chattanooga is a brand of DJO Global Inc. DJO Global provides solutions for musculoskeletal health, vascular health and pain management. Our products help prevent injuries or rehabilitate after surgery, injury or from degenerative disease, enabling patients to regain or maintain their natural motion. Visit www.DJOGlobal.eu



AIRCAST

chattanooga

CMR

Compex

djosurgical

DONJOY

Dr. Comfort


Empi

PROCARE

The Physical Therapy Solution for the Veterinary Market

Chattanooga, the world's largest manufacturer of rehabilitation equipment for treating musculoskeletal, neurological and soft tissue disorders, proudly brings their innovative products to the animal market.

From the original Hydrocollator HotPacs and heating units that Chattanooga Group pioneered 60 years ago, to the revolutionary advancements found in the Intellect Vet system, no other company knows physical therapy and rehabilitation like Chattanooga.



Clinical Indications for Dogs and Horses

Electrotherapy

- Muscle strengthening
- Wound healing
- Edema
- Pain control
- Increase circulation
- Fracture healing

Ultrasound Therapy

- Increase Circulation
- Wound healing
- Increase soft tissue mobility
- Decrease tendon and ligament inflammation

Laser Therapy

- Increase circulation
- Wound healing
- Tendonitis
- Pain control
- Decrease tendon and ligament inflammation



ORDERING INFORMATION

Part Number Description

Intellect Vet

2756 Intellect Vet Color Combo

Standard Accessories

27465 Patient Data Card (1x)

27335 5 cm² Ultrasound Applicator42192 Dura-Stick[®] Plus 5 cm round electrodes (4x)

4248 Ultrasound gel

28133 User Manual on CD

Optional Accessories

2767 Battery Module

2766 Laser Therapy Module

2770 Channel 3/4 Stimulation Module

27333 1 cm² Applicator27334 2 cm² Applicator27336 10 cm² Applicator

27133 Intellect Advanced Carrying Bag

2775ASY Therapy System Cart

79977 HiVolt Probe

57007 Microcurrent Probe

27525 Laser Protection Glasses (2x)

27900K Remote Interlock Kit

10648 Nylatex Wrap

79967 6 x 8 cm carbon electrodes (4x)

79970 6 x 8 cm sponges (4x)

INTELECT[®] VET

Intellect Vet is a revolutionary product for veterinarians who provide physical therapy to animals by offering four therapy modalities in one system:

- Electrical stimulation
- Ultrasound
- Combination stim and ultrasound
- Laser therapy

Intellect Vet is designed to be fully functional in the clinic or to be used in an environment where mobility is a must. An optional battery module is available, allowing therapy to be performed away from outlets. With the optional carry bag, transport and storage of the unit and accessories are easy.

FEATURES

- 5 clinical waveforms
- 2 independent electrotherapy channels
- High resolution color display of all screens
- Over 100 Clinical Protocols
- Over 100 User-defined protocol slots
- Impressive clinical library which includes full color graphic anatomical and pathological libraries on vibrant LCD screen
- 1 & 3 Mhz ultrasound frequencies
- Pulsed & Continuous Modes
- Ergonomically designed sound heads in 1 cm², 2 cm², 5 cm² & 10 cm² sizes
- Sound head coupling features with both visual & audio feedback
- 3 year warranty










TECHNICAL SPECIFICATIONS

Mains Power:	100-240 VAC, 50-60 Hz
Weight:	3.2 kg
Dimensions*:	33 x 24 x 29 cm
Electrical Safety Class:	Class I, Electrotherapy Type BF, Ultrasound Type B
Safety Tests:	IEC 60601-1 IEC 60601-1-2 IEC 60601-2-5 IEC 60601-2-10

* Dimensions are listed as L x W x H

WAVEFORMS

	INTERFERENTIAL <ul style="list-style-type: none"> • Four Pole • Manual Vector Scan • Two Pole (Premodulated)
	RUSSIAN <ul style="list-style-type: none"> • Continuous • Cycle Timed
IONTOPHORESIS	
	VMS <ul style="list-style-type: none"> • Continuous • Cycle timed • Burst
	HIGH VOLTAGE PULSED CURRENT
	ULTRASOUND <ul style="list-style-type: none"> • 1 and 3 MHz



ORDERING INFORMATION

Part Number Description

Laser Module

2766 Laser Therapy Module

Standard Accessories

27525 Laser Protection Glasses (2x)

28133 User Manual on CD

Optional Accessories

27900K Remote Interlock Kit

Laser Applicators

Singles

27799 670nm LED Diode 10mW

27803 850nm Laser Diode 40mW

27840 850nm Laser Diode 100mW

27804 850nm Laser Diode 150mW

27841 850nm Laser Diode 200mW

27805 820nm Laser Diode 300mW

9 Diode Cluster

27810 9 Diode Cluster Applicator, 290mW Total: 5 x 850nm 50mW Lasers, 4 x 670nm 10mW LED

27811 9 Diode Cluster Applicator, 540mW Total: 5 x 850nm 100mW Lasers, 4 x 670nm 10mW LED

27812 9 Diode Cluster Applicator, 1040mW Total: 5 x 850nm 200mW Lasers, 4 x 670nm 10mW LED

13 Diode Cluster

27813 13 Diode Cluster Applicator, 265mW Total: 3 x 850nm 50mW Lasers, 7 x 670nm 10mW LED, 3 x 950nm 15mW LED

27814 13 Diode Cluster Applicator, 415mW Total: 3 x 850nm 100mW Lasers, 7 x 670nm 10mW LED, 3 x 950nm 15mW LED

27816 13 Diode Cluster Applicator, 715mW Total: 3 x 850nm 200mW Lasers, 7 x 670nm 10mW LED, 3 x 950nm 15mW LED

19 Diode Cluster

27815 19 Diode Cluster Applicator, 325mW Total: LED only, 6 x 670nm 10mW LED, 7 x 880nm 25mW LED, 6 x 950nm 15mW LED

33 Diode Cluster

27809 33 Diode Cluster Applicator, 565mW Total: LED only, 12 x 670nm 10mW LED, 13 x 880nm 25mW LED, 8 x 950nm 15mW LED

27802 33 Diode Cluster Applicator, 690mW Total: 5 x 850nm 50mW Lasers, 12 x 670nm 10mW LED, 8 x 880nm 25mW LED, 8 x 950nm 15mW LED

27807 33 Diode Cluster Applicator, 940mW Total: 5 x 850nm 100mW Lasers, 12 x 670nm 10mW LED, 8 x 880nm 25mW LED, 8 x 950nm 15mW LED

27808 33 Diode Cluster Applicator, 1440mW Total: 5 x 850nm 200mW Lasers, 12 x 670nm 10mW LED, 8 x 880nm 25mW LED, 8 x 950nm 15mW LED

LASER MODULE




FEATURES

Low Level Laser Therapy

- Display of energy in choice of joules or joules/cm²
- Continuous and pulsed treatment modes with fully variable pulsed frequencies
- Real time feedback of energy delivered
- Wide variety of LED, Laser Diodes, & Cluster Probes available

The laser module has to be connected to the Intellect Vet color combo

 Visible and invisible laser radiation. Avoid direct exposure to beam. Class 3B Laser Product.

TECHNICAL SPECIFICATIONS

Weight: 1 kg (laser module only)

Dimensions*: 26 x 21 x 4 cm

Electrical Safety Class: Class 1, Type B

Safety Tests: IEC 60601-1, IEC 60825-1, EN 60601-2-22

* Dimensions are listed as L x W x H



BATTERY AND CHANNEL 3/4 STIM MODULES



Battery Module



Channel 3/4 Stimulation Module



FEATURES

Battery Module

- Nickel Metal Hydride battery self contained in packs of 10 cells
- Battery automatically and continually charges when Intellect Vet unit is plugged into mains supply

Channel 3/4 Stimulation Module

- 4 independent channels available on both Intellect Advanced Electrotherapy and Combination devices to upgrade unit from 2 channels

ORDERING INFORMATION

Part Number Description

Battery Module

2767 Battery Module

Standard Accessories

28133 User Manual on CD

4 Channel Stim Module

2770 Channel 3/4 Stimulation

27019 Dual Connector Cover

28133 User Manual on CD

Optional Accessories

79967 6 x 8 cm Carbon Electrodes (4x)

79970 6 x 8 cm Sponges (4x)

10648 Nylatex Wrap

TECHNICAL SPECIFICATIONS

Channel 3/4 Stimulation

Weight: 0.5 kg

Dimensions*: 30 x 21 x 4 cm

Electrical Safety Class: Class I, Type BF

Safety Tests: IEC 60601-1, IEC 60601-1-2, IEC 60601-2-10

Battery Module

Weight: 1.3 kg

Dimensions*: 30 x 21 x 4 cm

Electrical Safety Class: Class I, Type BF

Safety Tests: IEC 60601-1, IEC 60601-1-2

* Dimensions are listed as L x W x H

